

Deluxewipe A+

ALCOHOL BASED SURFACE DISINFECTANT WIPES



- Only 60 sec contact time
- Wide effective spectrum
- Powerful action
- Lemon aroma



High quality wet wipes designed for the disinfecting of alcohol resistant clinical surfaces.

The powerful alcohol based formula is effective against difficult to kill micro-organisms including Norovirus and Tuberculosis.

For the surface disinfection of non-invasive medical devices. Specially formulated for use in the medical and dental environment.

Busy dental practices need fast, effective solutions for maintaining hygiene between patients. Deluxewipes A+ from ViroSafe combine a high powered alcohol-based disinfectant solution with the convenience of pre-saturated disinfectant wipes. Designed for rapid cleaning and disinfection, Deluxewipes A+ deliver powerful antimicrobial performance in a simple, ready-to-use format.

Deluxewipes A+ provide broad-spectrum disinfection in just 60 seconds, ensuring reliable protection against bacteria, yeast and enveloped viruses. Their formulation is also effective against harder-to-kill microorganisms such as Mycobacteria, Norovirus, and Tuberculosis, giving dental professionals complete confidence in their infection control routines.

Each wipe is saturated with a 63% alcohol solution that ensures rapid evaporation for a streak-free finish. Deluxewipes A+ leave surfaces hygienically ready for use almost instantly - minimising downtime in busy clinical workflows.

Specially developed for use on alcohol-resistant materials, Deluxewipes A+ are perfectly suited for hard wearing surfaces commonly found in dental environments including worktops and equipment casings. They offer reliable disinfection without compromising the integrity of surfaces.

Enhanced with a subtle lemon fragrance, Deluxewipes A+ create a fresh, pleasant environment free from the sharp chemical odours often associated with alcohol disinfectants.

Packs of 160 wipes (175 x 155mm each). Refillable and reusable dispenser tubs for easy access and storage. Economical refill packs to reduce waste, costs and overall plastic consumption.

Instructions for Use.

Place refill bag in dispenser tub. Feed the centre wipe through the hole in the cap and replace. Using sufficient wipes for each application wipe the clean surfaces to be disinfected so that they are completely moistened and let the disinfectant take action over 60 seconds. Allow surface moisture to evaporate. For professional use only.

Composition.

Ethanol (64-17-5) 63% w/w (630g/kg), Aqua, Perfumes (Contains Linalool, Citral, Geraniol, Alpha-Isomethyl ionone).

Storage.

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Store in a well-ventilated place. Keep container tightly closed. Keep cool.

Contraindications.

Alcohol is not recommended for use on dental/medical instruments. Do not use on surfaces sensitive to alcohol.

Disposal.

Dispose of contents/container in accordance with local regulations. Do not discharge into drains or rivers. Do not flush. Non maceratable.

Efficacy and test standards



MICROORGANISM	TIME
Bactericidal	60 sec
EN13727	
EN16615	
Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, MRSA, VRE	60 sec
Yeasticidal	60 sec
EN13624	
EN16615	
Candida albicans	60 sec
Enveloped viruses	60 sec
EN14476	
Hepatitis B virus (HBV), Hepatitis C virus (HCV), HIV, Measles virus, Hepatitis Delta virus (HDV), Influenza virus, Coronavirus, Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Poxviridae	60 sec
Non enveloped viruses	60 sec
EN14476	
Adenovirus 5, Murine norovirus	60 sec
Mycobactericidal	60 sec
EN14348	
Mycobacterium avium	60 sec
Tuberculocidal	60 sec
EN14348	
Mycobacterium terrae	60 sec

Available in Dispenser tubs and refill packs 160 wipes 175 x 155mm

VIRS3007	Deluxewipes A+ Tub
VIRS3008	Deluxewipes A+ Refill
VIRS3009	Kit (6 Refills + 1 Empty Tub)
VIRS3010	Empty Tub

UFI: EE00-VOC9-F00M-WVDU



Manufactured in the United Kingdom by:
Topdental (Products) Ltd. Unit 1, Holmfield Industrial Estate, Halifax, W. Yorks, HX2 9TN. UK
Contact: info@virosafe.com, www.virosafeuk.com

* Always test an inconspicuous area first

v1.0 2508