

Surfacesafe F+

ALCOHOL FREE SURFACE DISINFECTANT

- Only 60 sec contact time
- Cleans and disinfects - 2 in 1
- Ideal to use on sensitive surfaces*
- Modern, safe formulation
- Apple aroma



Ready to use spray for the cleaning and disinfecting of clinical surfaces*

Perfect to use on a wide range of different surface types including those sensitive to alcohol, like the dental chair.

For the surface disinfection of medical devices, including but not limited to, dental chairs, examination lights and patient transport. Do not use on invasive instruments. Specially formulated for use in the medical and dental environment.

Surfacesafe F+ from ViroSafe is a next-generation, alcohol-free surface disinfectant spray designed specifically for the needs of modern dental clinics. With rapid action, material safety, and a subtle, refreshing fragrance, Surfacesafe F+ delivers professional-grade disinfection you can trust every day.

Time is critical in a busy clinical environment. Surfacesafe F+ delivers powerful disinfection in just 60 seconds, ensuring rapid turnover between patients without compromising infection control. Proven effective against a broad spectrum of bacteria, Yeast and all enveloped viruses such as Coronavirus. Surfacesafe F+ provides dependable protection against cross-contamination in high-contact areas.

Unlike harsh alcohol-based sprays that can damage sensitive equipment, Surfacesafe F+ is almost universally compatible with a wide variety of materials. From dental chairs, plastics, and synthetic leather to stainless steel and glass, its gentle formulation protects your investment while ensuring surfaces remain spotless and hygienically safe*

Surfacesafe F+ contains a plant-derived detergent that enhances its cleaning performance while supporting a more sustainable approach to infection control. Tough on dirt and organic residues yet gentle on surfaces, this eco-conscious ingredient ensures cleanliness and safety go hand in hand.

Formulated with a light apple fragrance, Surfacesafe F+ leaves the treatment area feeling clean and refreshed - without the strong chemical odours often linked to disinfectants.

Instructions for Use.

Spray clean surfaces and objects to be disinfected so that they are completely moistened and wipe with a clean cloth to distribute. Let the disinfectant take action over 60 seconds then allow to air dry or wipe with clean tissue. Contains detergent. For professional use only.

Composition.

Contains Aqua, Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-16)) (68424-85-1) 0.48% w/w (4.8g/kg), Didecyldimethylammonium Chloride (7173-51-5) 0.22% w/w (2.2g/kg), Perfumes (contain Amyl cinnamal, Limonene, Alpha-Isomethyl ionone).

Handling and Storage.

Ensure there is sufficient ventilation of the area. The wearing of protective gloves is recommended. Store in a well-ventilated place. Keep container tightly closed. Keep cool.

Contraindications.

Mild steel, aluminium, copper.

Disposal.

Dispose of contents/container in accordance with local regulations.
Do not discharge liquid into drains or rivers.

Harmful to aquatic life with long lasting effects. Avoid release to the environment. IF SWALLOWED: Immediately call a doctor. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.



Efficacy and test standards

MICROORGANISM	CONTACT TIME
Bactericidal	60 sec
EN13727	
EN16615	
Pseudomonas aeruginosa	60 sec
Staphylococcus aureus	60 sec
Enterococcus hirae	60 sec
Yeast	60 sec
EN13624	
EN16615	
Candida albicans	60 sec
Enveloped viruses	60 sec
EN14476	
EN16777	
Vaccinia virus	60 sec
Hepatitis B virus (HBV)	60 sec
Hepatitis C virus (HCV)	60 sec
Hepatitis Delta virus (HDV)	60 sec
HIV	60 sec
Influenza virus	60 sec
Coronavirus	60 sec
Herpesviridae	60 sec
Filoviridae (Ebola)	60 sec
Rubella virus	60 sec
Rabies virus	60 sec
Measles virus	60 sec
Poxviridae	60 sec

REF Available in 1 Litre, 2.5 Litre, 5 Litre

VIRS4001 Surfacesafe F+ 1 litre
VIRS4002 Surfacesafe F+ 2.5 litre
VIRS4003 Surfacesafe F+ 5 litre

UFI: 2520-0057-K00H-5Q24



Manufactured in the United Kingdom by:
Topdental (Products) Ltd. Unit 1, Holmfield
Industrial Estate, Halifax, W. Yorks, HX2
9TN, UK
Contact: info@viroSAFE.com, www.viroSAFEuk.com

* Always test an inconspicuous area first

v1.0 2508